

MW45



Nasopharyngeal Swab Collection Simulator

For training to healthcare professionals who conduct virus testing for the first time



FEATURES

- 1 | When swab reaches the right position in the nasal cavity, the swab will be pigmented in blue
- 2 | Ways to use tongue depressor with swabs can be practiced.
- 3 | The shape of the left and right nasal cavity is different.
- 4 | The product comes with a cross sectional model.

SKILLS

- | Training to healthcare professionals who conduct virus testing for the first time.

DESCRIPTIONS

SPECIFICATIONS

Size: W42 x D23 x H38 cm/W16.8 x D9.2 x H15.2 in Weight: 2.15 kg / 4.73 lbs

MATERIALS

Soft resin, Hard resin
Latex free

SET INCLUDES

- 1 main model
- 1 cross sectional model
- 1 syringe (10ml)
- 1 disposable cup
- 10 cotton swab blue (set of 10)
- 1 instruction manual
- 1 storage case
- 6 replacement pads

REPLACEMENT PARTS

- 11426-010 replacement pad (set of 5)
- 11402-010 cotton swab blue (set of 10)

MW60

Tracheostomy Simulator

Simulator for percutaneous tracheostomy, puncture and incision of cricothyroid ligament



FEATURES

- 1 | Anatomically correct cricoid cartilage and thyroid cartilage with life-like senses of touch.
- 2 | Multiple use possible by rotation of neck skins.
- 3 | Adhesive discs to hold the manikin head in place during the procedure.

SKILLS

- | Percutaneous tracheostomy
- | Puncture and incision of cricothyroid ligament

DESCRIPTIONS

SPECIFICATIONS

Size: W13.5 x D24 x H12.5cm Weight: 1.7kg
W17.7 x 7.9 x H9.8 inch

SET INCLUDES

- 1 main body
- 1 skin for tracheostomy(set of 3)
- 1 subcutaneous tissue for
- 1 tracheostomy(set of 3)
- 1 cricoid cartilage
- 1 instruction manual

MATERIALS

soft resin, Hard resin
Latex free

REPLACEMENT PARTS

- 11442-010 Skin for tracheostomy (set of 3)
- 11442-020 subcutaneous tissue for tracheostomy (set of 3)