

PH-26 | 41503-000

## ORINS Thyroid Phantom ITS



A phantom by the ORINS standards



SHOW MORE!



### FEATURES

- | Oak Ridge Institute for Nuclear Studies type phantom for measurement of thyroid radionuclide uptake
- | Cavities for iodine-131 are prepared in the neck phantom

### APPLICATIONS

- | SPECT

### DESCRIPTIONS

#### SET INCLUDES

- 1 petroleum jelly
- 1 screwdriver
- 1 storage case
- 1 manual

#### MATERIALS

Acrylic resin

#### SPECIFICATIONS

Phantom size:  
12.5 dia. x 12.5 (H) cm/4.9dia. x 4.9(H) in

PH-29 | 41540-030

## ECT Hot Cold Phantom SP-6



Volumetric measurement phantom for PET/SPECT



SHOW MORE!



### FEATURES

- | Five sphere containers with different sizes can be filled with RI solution
- | Volume of sphere phantoms are:  
50 mm/2 in (100%), 80%, 60%, 40% and 20%

### DESCRIPTIONS

#### SET INCLUDES

- 1 phantom
- 1 storage case

#### MATERIALS

Acrylic resin

#### SPECIFICATIONS

Phantom size:  
21 dia. x 16 (H) cm / 8.2 dia. x 6.2 (H) in

### APPLICATIONS

- | SPECT, PET