


## MW82 Femoral Arterial Puncture Trainer

Hands-on femoral arterial access training both with ultrasound-guided and landmark-guided procedures.

### Features

- Realistic arterial pulse
- Palpable ASIS and pubic bone landmarks
- Simulated reverse blood flow
- No visible puncture marks at the injection site
- Leak-proof connections



### Skills

- Ultrasound-guided femoral artery puncture
- Landmark-guided femoral artery puncture
- Blood sampling

# Femoral Arterial Puncture Trainer

MW82

11486-000

## Set Includes

groin manikin	1
circulation pump (AC100V-240V)	1
artery tube	2
vein tube	2
cotton swab red	10
instruction manual	

## Specifications

	Groin manikin
dimensions	W33 × D35 × H15 cm (manikin) W20 × D16 × H13 cm (pump)
weight	3.4 kg (manikin) 1.5kg (pump)
materials	Soft resin, Latex free



## Recommended devices

Needle size: 23G

Replacement parts



replacement skin  
11486-030



artery tube (a set of 4)  
11486-010



vein tube (a set of 4)  
11486-020



cotton swab red (set of 10)  
11388-400

Related products




An essential simulator to increase confidence and competence

Arterial Puncture Training Wrist  
M99 11351-000



Master Ultrasound-Guided Radial Arterial Puncture

Ultrasound Guided Arterial Puncture Training Pads  
11351-110



A sequence of uninterrupted simulation experience for specialist training

CVC Insertion Simulator II & III  
M93UB, U, C 11347-300, 500, 510

**KYOTO KAGAKU**

<https://www.kyotokagaku.com/>  
[rw-kyoto@kyotokagaku.co.jp](mailto:rw-kyoto@kyotokagaku.co.jp)

2026.03

### Kyoto Head Office / Factory / Show Room

15 Kitaneokoya-cho Fushimi-ku Kyoto, Japan 612-8388  
Telephone:+81-75-605-2510 Fax:+81-75-605-2519

### Asia and other regions

#### Overseas Div. / MEM Training Center

Second fl. NREG Hongo-3-chome bldg., 3-26-6 Hongo,  
Bunkyo-ku, Tokyo 113-0033, Japan  
Telephone:+81-3-3817-8071 Fax:+81-3-3817-8075

### All American regions

#### Kyoto Kagaku America Inc. / Show Room

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA  
Telephone:+1-310-325-8860 Fax:+1-310-325-8867  
Email : info@kkamerica-inc.com

### Europe, Russia & Africa

#### Kyoto Kagaku Europe GmbH / Training Room

De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany  
Telephone:+49-69-5060-28160