

US-23

41953-000




NEW

Rectal and Bladder Ultrasound Phantom

Ultrasound training cubes for assessment of residual urine and excretion

Product Supervision

Hiromi Sanada, RN, PhD, WOCN, FAAN President,
Ishikawa Prefectural Nursing University
Masaru Matsumoto, Asst. Professor
Ishikawa Prefectural Nursing University





FEATURES

- Variation of residue excretion in the rectum for ultrasound assessment:
hard, soft, combination and empty
- Three types of urinal bladder for assessment and urine volume measurement:
small, middle, large with a balloon catheter



SKILLS

- Ultrasound assessment of residue excretion in the rectum
- Measurement of urine volume
- Ultrasound assessment of the urinal bladder and the prostate
- Assessment of indwelling bladder catheter



A: Small bladder / Hard excretion



B: Middle size bladder / Soft excretion



C: Large bladder with a balloon catheter / Hard and soft excretion / empty rectum



DESCRIPTIONS

SPECIFICATIONS

Size: W15 x D15 x 20H cm each
W5.9 x D5.9 x H3.4 in. each

Weight: 2.0kg each
0.9lbs each

SET INCLUDES

- 3 types of ultrasound phantom cube
- 1 talcum powder
- 1 storage case
- 1 instruction manual

* Ultrasound scanner is not included.